

# HotStart Taq DNA Polymerase

*For research use only*

**Geneaid**



CERTIFICATE NO. QAIC/TW50077  
**ISO 9001:2008 QMS**

**Catalogue Numbers:** HS005, HS050, HS250, HSN005, HSN050, HSN250

**Quantity:** 50 U, 500 U, 2500 U

**Unit Definition:** 1 U incorporates 10 nmole of dNTP into acid insoluble products in 30 minutes at 74°C

**Storage Buffer:** 20 mM Tris-HCl pH8.0, 0.1 mM EDTA, 1 mM DTT, 1.0% Triton X-100, 50% Glycerol

**10X PCR Buffer II:** 200 mM Tris-HCl, pH 8.1, 40 mM MgCl<sub>2</sub>, 100 mM KCl, 1.0% Triton X-100

**MgCl<sub>2</sub>:** premixed with 10X PCR Buffer

**Storage:** -20°C for extended periods

**3' to 5' Exonuclease Proofreading Ability:** NO

**5' to 3' Exonuclease Activity:** YES

## Introduction

HotStart Taq DNA Polymerase is a thermostable, inhibitor modified form of recombinant Taq DNA Polymerase with ultra low DNA content. It is inactive at temperatures below 45°C, but is activated at 95°C (1 minute). HotStart Taq DNA Polymerase is ideal for "Hot Start" qPCR, because the enzyme remains inactive during qPCR set-up. Since it is inactive at low temperatures, it does not elongate non-specific primer-template hybrids that may form at low temperatures. The enzyme will amplify DNA targets with high specificity, sensitivity, and yield.

## Quality Control

HotStart Taq DNA Polymerase is tested on a lot-to-lot basis. Quality control assays include activity test, PCR, endonuclease activity, exonuclease activity.

## Caution

During the procedure, always wear a lab coat, disposable gloves, and protective goggles

Component	HS005/HSN005	HS050/HSN050	HS250/HSN250
HS Taq DNA Polymerase	10 µl (50U)	100 µl (500U)	100 µl x 5 (2500U)
10X PCR Buffer II	200 µl	2 ml	2 ml x 5
dNTPs (10 mM each)	NA/20 µl	NA/200 µl	NA/1 ml

## Advantages

**Simplified "Hot Start" Technique:** The modified Taq DNA Polymerase is inactive at room temperature but is readily activated during the first high temperature step of PCR.

**Fast Activation (95°C for 1 minute):** Other Hot Start Taq requires 4 to 10 minutes of heating to activate the enzyme. Extended heating time reduces the duration of Taq activity and can delay the cycle threshold (Ct).

**Increased PCR Specificity:** The enzyme will not amplify non-specific, low-temperature primer-template hybrids.

**Increased PCR Sensitivity & Yield:** The enzyme leads to a more specific product.

**Improved Performance of Multiplex PCR:** Eliminates non-specific priming, resulting in higher specificity/sensitivity, and lower non-specific background signals.

**Ideal for High Throughput or Automated PCR:** Reactions can be set up at room temperature.

**Reduces Risk of Contaminant Carry-over:** Allows "Hot Start" without additional pipetting/handling steps or additional reaction components such as wax.

## HotStart Taq PCR Reaction Set Up


All reaction components should be prepared on ice or in a biocooler then quickly transferred to a thermal cycler pre-heated to 95°C. For added convenience and to prevent pipetting errors, combine the reaction components in a PCR master mix then aliquot several reactions in parallel PCR tubes. Storage Buffer contains 50% Glycerol so HotStart Taq DNA Polymerase does not require thawing when removed from -20°C storage.

1. Allow components to thaw then vortex gently and spin down.
2. Remove HotStart Taq DNA Polymerase from -20°C storage and spin down.
3. Add the following components to a PCR tube or tubes.

Reaction Components	Volume (µl)	Amount
Template DNA	Variable	<500 ng
Upstream Primer (10 µM)	0.2-1.0 µl	0.1-0.5 µM
Downstream Primer (10 µM)	0.2-1.0 µl	0.1-0.5 µM
dNTP Mix (10 mM each)	0.4 µl	200 µM
10X PCR Buffer II	2 µl	4 mM MgCl <sub>2</sub>
HS Taq DNA polymerase	0.2 µl	1 Unit
ddH <sub>2</sub> O	Add to 20 µl	

## General Thermal Cycling Program

Please read the entire instruction manual prior to starting the thermal cycling program.

Initial Denaturation	<ul style="list-style-type: none"><li>• 95°C.</li><li>• 2-5 minutes.</li></ul> <p>*Denaturation condition varies depending on thermal cycler and tube.</p>	
Step 2 Denaturation	<ul style="list-style-type: none"><li>• 95°C.</li><li>• 20-40 seconds.</li></ul>	
Step 3 Annealing	<ul style="list-style-type: none"><li>• 55°C-60°C.</li><li>• 30-60 seconds.</li><li>• Optimal annealing temperature is dependent on the T<sub>m</sub> value of the primers and the reaction condition.</li></ul>	
Step 4 Elongation	<ul style="list-style-type: none"><li>• 72°C</li><li>• 20 seconds (for amplicon &lt;500 bp).</li><li>• For large amplicon PCR, add 30 seconds for each kbp.</li></ul>	
Step 5 Final Elongation	<ul style="list-style-type: none"><li>• 72°C.</li><li>• 7 minutes.</li></ul>	
Step 6 Cooling	<ul style="list-style-type: none"><li>• 4°C.</li></ul>	

## Patent Disclaimer

Some applications in which this product can be used may be covered by patents issued and applicable in certain countries. Because purchase of this product does not include a license to perform any patented application, users of this product may be required to obtain a patent license depending on the particular application and country in which the product is being used.