

# ZEROprep™ DNA Reagent

For research use only

Store at room temperature

Catalogue Numbers

ZDR005  
ZDR050  
ZDR100

Quantity

500 µl  
50 ml  
100 ml

Geneaid



CERTIFICATE NO. QAIC/TW/50077  
ISO 9001:2008 QMS

## Introduction

ZEROprep™ DNA Reagent is designed for efficient release of DNA for direct use in PCR reactions without purification. A wide variety of samples are effectively homogenized in the reagent without any pre-treatment or subsequent bind, wash or elution steps. Simply place the sample in the reagent, follow the 2 step protocol and transfer the lysate to a PCR mix.

## Quality Control

ZEROprep™ DNA Reagent is tested on a lot-to-lot basis according to Geneaid's ISO-certified quality management system. DNA from a 1 mg tissue sample is lysed in ZEROprep™ DNA Reagent. A 5 µl aliquot of lysate is added directly into a 50 µl PCR mix.

## Advantages

- Use DNA directly in PCR reactions
- DNA purification is not required
- 15 minute 2 step protocol
- Wide variety of sample types (tissue, blood, plant, bacteria, yeast/fungus, virus)

## Applications

Direct use of DNA in PCR reactions, multiplex PCR, Real-time PCR

## Caution

During operation, always wear a lab coat, disposable gloves, protective goggles or (anti-fog) procedure mask.

## ZEROprep™ DNA Reagent Protocol Procedure

Please read the entire instruction manual prior to starting the Protocol Procedure

Sample	Procedure
Tissue	<ol style="list-style-type: none"><li>1. <b>Transfer 50 µl of ZEROprep™ DNA Reagent and 1 mg of tissue to a 1.5 ml microcentrifuge tube.</b></li><li>2. Incubate for 15 minutes at room temperature or 5-15 minutes at 80°C.</li><li>3. Mix by vortex then transfer a <b>2-5 µl aliquot</b> to a <b>20-50 µl PCR mix.</b></li></ol>
Plant Tissue <sup>1</sup>	<ol style="list-style-type: none"><li>1. <b>Transfer 200 µl of ZEROprep™ DNA Reagent and 5-25 mg of tissue to a 1.5 ml microcentrifuge tube.</b></li><li>2. Incubate for 15 minutes at room temperature or 5-15 minutes at 80°C.</li><li>3. Mix by vortex then transfer a <b>2-5 µl aliquot</b> to a <b>20-50 µl PCR mix.</b></li></ol>
Whole Blood, plasma, serum	<ol style="list-style-type: none"><li>1. <b>Transfer 100 µl of ZEROprep™ DNA Reagent and 5-10 µl of fluid sample to a 1.5 ml microcentrifuge tube.</b></li><li>2. Incubate for 15 minutes at room temperature.</li><li>3. Mix by vortex then transfer a <b>2-5 µl aliquot</b> to a <b>20-50 µl PCR mix.</b></li></ol>
Saliva	<ol style="list-style-type: none"><li>1. <b>Transfer 100 µl of ZEROprep™ DNA Reagent and 10 µl of saliva to a 1.5 ml microcentrifuge tube.</b></li><li>2. Incubate for 15 minutes at room temperature or 10 minutes at 80°C.</li><li>3. Mix by vortex then transfer a <b>2-5 µl aliquot</b> to a <b>20-50 µl PCR mix.</b></li></ol>
Bacteria <sup>2</sup>	<ol style="list-style-type: none"><li>1. <b>Transfer 100 µl of ZEROprep™ DNA Reagent and 1-5 µl of bacteria culture to a 1.5 ml microcentrifuge tube.</b></li><li>2. Incubate for 15 minutes at room temperature or 10-15 minutes at 80-90°C.</li><li>3. Mix by vortex then transfer a <b>2-5 µl aliquot</b> to a <b>20-50 µl PCR mix.</b></li></ol>

NOTE:  
<sup>1</sup>For plant species with high levels of polysaccharide inhibitors, increase the tissue amount by 2-3 times per volume of ZEROprep™ DNA Reagent. Plant tissue homogenization using a bead beating instrument or pestle and mortar with liquid nitrogen will facilitate DNA release.  
<sup>2</sup>*E. coli* can be efficiently lysed in ZEROprep™ DNA Reagent for 15 minutes at room temperature. However, to efficiently disrupt the bacteria cell wall of gram (+) bacteria, 3 hour incubation at room temperature or 10-15 minutes at 80°C is required.