

# Magnetic Beads gDNA Kit (Blood)

For research use only

## Catalogue Numbers

MB048

MB096

## Quantity

48 rxns

96 rxns

Geneaid



CERTIFICATE NO. QAIC/TW/50077  
ISO 9001:2008 QMS

## Introduction

The Magnetic Beads Genomic DNA Extraction Kit Blood was designed specifically for efficient genomic DNA purification from whole blood and buffy coat. DNA is bound to the surface of the magnetic beads and released using a proprietary buffer system. The Magnetic Beads Genomic DNA Extraction Kit Blood can be easily adapted to automated magnetic bead separation instruments and workstations. The purified DNA can be used in qPCR and a variety of other downstream applications.

## Quality Control

The quality of the The Magnetic Beads Genomic DNA Extraction Kit Blood is tested on a lot-to-lot basis according to Geneaid's ISO-certified quality management system by isolating genomic DNA from a 300 µl whole blood sample.

## Advantages

- High Yield: 10 µg of Genomic DNA
- High Quality DNA: A260/A280 = 1.8-2.0
- Easily adapted to automated magnetic bead separation instruments and workstations
- Sample: 100 to 400 µl of whole blood and buffy coat (5 x 10<sup>6</sup> WBC)
- Operation time: within 20 minutes (manual)
- Storage: dry at room temperature (15-25°C) for up to 1 year, Protease should be stored dry at 2-8°C for up to 6 months

## Caution

During operation, always wear a lab coat, disposable gloves, protective goggles and (anti-fog) procedure mask.

## Components and Storage

Item	Volume	Product	Shipping	Storage
MB1 Buffer	2 ml	MB004	room temperature	dry at room temperature (15-25°C)
	30 ml	MB048		
	60 ml	MB096		
MB2 Buffer <sup>1</sup> (Add Isopropanol)	0.8 ml (1.1 ml)	MB004	room temperature	dry at room temperature (15-25°C)
	11 ml (14 ml)	MB048		
	22 ml (28 ml)	MB096		
Protease <sup>2</sup>	100 µl	MB004	room temperature	dry at 2-8°C
	1 ml	MB048		
	1 ml x 2	MB096		
MW1 Buffer	2 ml x 2	MB004	room temperature	dry at room temperature (15-25°C)
	45 ml	MB048		
	60 ml	MB096		
MW2 Buffer <sup>3</sup> (Add Ethanol)	1 ml (4 ml)	MB004	room temperature	dry at room temperature (15-25°C)
	12.5 ml (50 ml)	MB048		
	25 ml (100 ml)	MB096		
MB Magnetic Beads	220 µl	MB004	room temperature	dry at room temperature (15-25°C)
	2.5 ml	MB048		
	5 ml	MB096		
Elution Buffer	1 ml	MB004	room temperature	dry at room temperature (15-25°C)
	12 ml	MB048		
	30 ml	MB096		

<sup>1</sup>Add Isopropanol (see the bottle label for volume) to MB2 Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid Isopropanol evaporation.

<sup>2</sup>Protease is shipped at room temperature and should be stored dry at 2-8°C for up to 6 months

<sup>3</sup>Add absolute ethanol (see the bottle label for volume) to MW2 Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

# Magnetic Beads Genomic DNA Extraction Kit Protocol Procedure

## IMPORTANT BEFORE USE:

1. Vortex magnetic beads to ensure they are in suspension prior to initial use.
2. Be sure and allow magnetic beads to disperse completely during the binding, wash and elution steps.
3. Add Isopropanol (see the bottle label for volume) to MB2 Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid Isopropanol evaporation.
4. Add absolute ethanol (see the bottle label for volume) to MW2 Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

**Additional requirements:** absolute ethanol, microcentrifuge tubes, magnetic separator, isopropanol

**1. Transfer 100-400 µl of blood (WBC count less than  $2 \times 10^4$  cells/µl) or 100-400 µl of buffy coat (less than  $5 \times 10^6$  WBC) into a 1.5 ml microcentrifuge tube. Add 400 µl of MB1 Buffer and 20 µl of Protease. Mix well and incubate at 65°C for 5 minutes. During incubation, invert the tube occasionally.**

**2. Add 450 µl of MB2 Buffer (make sure isopropanol was added) to the sample and mix well by vortex.** Vortex the MB Magnetic Beads for 10 seconds prior to use to ensure the MB Magnetic Beads are in suspension. **Add 50 µl of MB Magnetic Beads.** Gently shake the tube for 5 minutes to ensure the MB Magnetic Beads disperse completely in the sample mixture. Place the tube in a magnetic separator for 30 seconds or until MB Magnetic Beads have pelleted. Remove and discard the supernatant.

**3. Add 600 µl of MW1 Buffer** and gently shake the tube for 1 minute. Place the tube in a magnetic separator for 30 seconds or until MB Magnetic Beads have pelleted. Remove and discard the supernatant. **Add 600 µl of MW2 Buffer (make sure ethanol was added)** and gently shake the tube for 1 minute. Place the tube in a magnetic separator for 30 seconds or until MB Magnetic Beads have pelleted. Remove and discard the supernatant. **Add 600 µl of MW2 Buffer (make sure ethanol was added)** and gently shake the tube for 1 minute. Place the tube in a magnetic separator for 30 seconds or until MB Magnetic Beads have pelleted. Remove and discard the supernatant.


**4. Incubate the tube at 65°C for 3 minutes to dry the MB Magnetic Beads. Add 50–200 µl of Elution Buffer.** Mix the sample by pipetting then incubate at room temperature for 3 minutes. During incubation, keep the MB Magnetic Beads in suspension by mixing. Place the tube in a magnetic separator for 30 seconds or until MB Magnetic Beads have pelleted. Carefully transfer the supernatant containing the purified DNA to a clean 1.5 ml microcentrifuge tube.

## Magnetic Beads Genomic DNA Extraction Kit Blood Functional Test Data



Sample	Volume	No. of WBC	Yield	260/280	260/230
Whole	100 µl	$5.0 \times 10^5$	2.5-3.0 µg	>1.8	>1.2
	200 µl	$1.0 \times 10^6$	5.5-6.0 µg	>1.8	>1.6
Blood	300 µl	$1.5 \times 10^6$	7.5-8.0 µg	>1.8	>1.7
	400 µl	$2.0 \times 10^6$	9.5-10.0 µg	>1.8	>1.8

**Figure 1.** 5 µl of extracted DNA product from a 100 µl eluate was analyzed on a 0.8% agarose gel.



**Figure 2.** Extracted DNA products were used as a DNA template for amplifying partial human ACTB gene. 3 µl of PCR product was loaded in each well. A 250 bp ACTB gene fragment was successfully amplified from each DNA product.

M1/M2: Geneaid 1 Kb DNA Ladder

Lane 1-3: Extracted DNA from 100 µl blood sample

Lane 4-6: Extracted DNA from 200 µl blood sample

Lane 7-9: Extracted DNA from 300 µl blood sample

Lane 10-12: Extracted DNA from 400 µl blood sample