

Presto™ Food DNA Extraction Kit Quick Protocol

For research use only

Catalogue Number

FGD004, FGD100, FGD300

Geneaid



Instruction Manual Download

When using this product for the first time, or if you are unfamiliar with the procedure, please scan the QR code and download the complete instruction manual.

Important Before Use

1. Add ddH₂O to Proteinase K (see the bottle label for volume) then vortex to ensure it is completely dissolved. Check the box on the bottle. For extended periods, the ddH₂O and Proteinase K mixture should be stored at 4°C. Use only fresh ddH₂O as ambient CO₂ can quickly cause acidification.
2. Add absolute ethanol (see the bottle label for volume) to Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

Additional Requirements

2.0 ml microcentrifuge tubes, 1.5 ml microcentrifuge tubes, standard vortex, 60°C incubator, absolute ethanol.

Presto™ Food DNA Extraction Kit Protocol

1. Sample Lysis

Homogenize **200 mg of food sample** using one of the following methods: **A.** Add liquid nitrogen to a mortar and grind the food samples thoroughly using a pestle. Transfer the food powder to a 2 ml centrifuge tube. **B.** Transfer the food samples to a 2 ml centrifuge tube containing ceramic beads or stainless steel beads then homogenize the sample with a TissueLyser, Disruptor Genie or similar. **C.** Transfer 200 µl of liquid samples directly to a 2 ml centrifuge tube.

Add 1 ml of **FD1 Buffer** and **10 µl of Proteinase K** to the microcentrifuge tube containing the homogenized food sample then vortex briefly.

NOTE: For samples which swell from liquid absorption (e.g., starches, cornflakes), add 2 ml of FD1 Buffer to ensure the sample is completely submerged.

Incubate the sample in a 60°C incubator with shaking for 30 minutes. If the incubator doesn't have a shaking function, manually invert the tube every 5 minutes.

NOTE: Preheat the required **Elution Buffer** (200 µl per sample) to 60°C for DNA elution.

2. PCR Inhibitor Removal

Cool the sample on ice or at 4°C to room temperature to facilitate inhibitor precipitation. Centrifuge the tubes at 2,500 x g for 5 minutes at room temperature. Transfer **600-700 µl of clear supernatant to a clean 1.5 ml microcentrifuge tube.**

NOTE: Depending on the type of food, the supernatant may be colored. If a semi-solid layer floats on top of the supernatant, pierce it with a pipette and transfer only the clear supernatant to a clean 1.5 ml microcentrifuge tube, being careful **NOT** to touch any precipitate/pellet from the bottom of the tube.

Add **500 µl of Chloroform** then vortex for 15 seconds. Centrifuge at 16,000 x g for 10 minutes at room temperature. If the supernatant is not clear, centrifuge again for 5 minutes. Carefully transfer **500 µl of the upper, aqueous phase** to a new 1.5 ml microcentrifuge tube.

3. DNA Binding

Add **500 µl of FD2 Buffer** and mix **IMMEDIATELY** by shaking the tube vigorously for 5 seconds.

NOTE: For extracting small DNA fragments (100-200 bp) from highly processed food samples, add 1 ml of FD2 Buffer to 250 µl of the sample instead.

Place a **GD Column in a 2 ml Collection Tube.** Transfer **700 µl of sample mixture** to the GD Column then centrifuge at 16,000 x g for 1 minute at room temperature. Discard the flow-through. Place the **GD Column back in the 2 ml Collection Tube.** Transfer all of the remaining sample mixture to the GD Column then centrifuge at 16,000 x g for 1 minute at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube.

4. Wash

Add **400 µl of W1 Buffer.** Centrifuge at 16,000 x g for 30 seconds at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube. Add **600 µl of Wash Buffer (make sure absolute ethanol was added).** Centrifuge at 16,000 x g for 30 seconds at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube. Centrifuge at 16,000 x g for 3 minutes at room temperature to dry the column matrix.

Optional Residual Pigment Removal Step

If pigments remain on the column, perform this optional step.

Following Wash Buffer addition, add 500 µl of absolute ethanol to the GD Column. Centrifuge at 16,000 x g for 30 seconds at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube. Centrifuge at 16,000 x g for 3 minutes at room temperature to dry the column matrix.

5. Elution

Transfer the dry **GD Column** to a new 1.5 ml microcentrifuge tube. Add **100 µl of preheated Elution Buffer**¹, TE² or water³ into the CENTER of the column matrix. Let stand for at least 2 minutes to allow Elution Buffer, TE or water to be completely absorbed. Centrifuge at 16,000 x g for 2 minutes at room temperature to elute the purified DNA.

¹If maximum DNA yield is required, use 200 µl of Elution Buffer (DNA concentration will be diluted). Ensure that Elution Buffer is added into the center of the GD Column matrix and is completely absorbed.

²Using TE (10 mM Tris-HCl, 1 mM EDTA, pH8.0) for elution is beneficial as EDTA preserves DNA for long term storage. However, EDTA will affect PCR and other sensitive downstream applications. Ensure that TE is added into the center of the GD Column matrix and is completely absorbed.

³If using water for elution, ensure the water pH is between 7.0 and 8.5. ddH₂O should be fresh as ambient CO₂ can quickly cause acidification. Ensure that water is added into the center of the GD Column matrix and is completely absorbed. DNA eluted in water should be stored at -20°C to avoid degradation.


Components

	FGD004	FGD100	FGD300
FD1 Buffer	10 ml	200 ml	200 ml x3
FD2 Buffer ²	2 ml x2	100 ml	100 ml x3
Proteinase K ¹ (Add ddH ₂ O)	1 mg (0.1 ml)	11 mg (1.1 ml)	11 mg x3 (1.1 ml)
W1 Buffer	2 ml	45 ml	130 ml
Wash Buffer ² (Add Ethanol)	1 ml (4 ml)	25 ml (100 ml)	50 ml (200 ml)
Elution Buffer	1 ml	30 ml	75 ml
GD Columns	4	100	300
2 ml Collection Tubes	4	100	300

¹ Add ddH₂O to Proteinase K (see the bottle label for volume) then vortex to ensure it is completely dissolved. Check the box on the bottle. Once it is dissolved completely, centrifuge for a few seconds to spin the mixture down. For extended periods, the ddH₂O and Proteinase K mixture should be stored at 4°C. Use only fresh ddH₂O as ambient CO₂ can quickly cause acidification.

² Add absolute ethanol (see the bottle label for volume) to Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

Presto™ Food DNA Extraction Kit Functional Test Data



Sample	ng/µl	260/280	260/230	Yield (µg)
1	95.70	1.89	2.31	09.6
2	110.04	1.90	2.30	11.0
3	108.50	1.89	2.29	10.9

Figure 1. Total DNA was purified from 200 mg of sausage using the Presto™ Food DNA Extraction Kit. 10 µl aliquots of extracted DNA from a 100 µl eluate was analyzed on a 0.8% agarose gel. M = Geneaid™ 1 Kb DNA Ladder